

	Technical Data Sheet	Pumpmatic Liquid Dispensing System
	10 ml Serological Pumpmatic	PN 110 and 210

Uses For tissue culture or any non-toxic solution for aspirating/dispensing

Materials: Syringe	Polypropylene; plunger non-latex rubber
Luer Adapter	Polystyrene
Pipet	Polystyrene
Individual Packaging	63 grams coated paper + 0.052mm PET/PE transparent film pouch
Carton Packaging	Corrugated Cardboard box
--Box Size - 110:	19 x 16 x 5 Inches
--Box Size - 210:	15 x 8 x 9 Inches

Regulatory:	FDA CGMP Practices	
	FDA Classification Name	Pipette, Diluting
	Product Code	GGY
	Device Class	1
	Regulation Number	864.6160
	Registered Establishment No	1525728

Related Catalog Numbers:		
	Sterile 110	Individually Wrapped, 25 per bag, 100 per case
	Non-Sterile 210	Bulk, 25 per bag, 100 per case

#### General Dimensions

Overall Length w/ Syringe	14.00 (inches)
Overall Length (pipet only)	10.00 (inches)
Pipet Size	10ml
Increments	0.1
Accuracy within $\pm 2\%$ at full volume	
Calibrated to deliver (TD-EX 20°C) their contents with blowout	

Sterility Method (110 only)	Electron Beam	13.2 - 28.0 kGy exposure
-----------------------------	---------------	--------------------------

